|  |   |    |
|--|---|----|
| MRT-Sicherheitsinformationen                             | 2 | DE |
| MRI Safety Information                                   | 3 | EN |
| Informations de sécurité IRM                             | 4 | FR |
| Informazioni sulla sicurezza RMI                         | 5 | IT |
| Información de seguridad sobre resonancia magnética (RM) | 6 | ES |
| Informações de segurança para RMN                        | 7 | PT |
| Информация по безопасности при МРТ                       | 8 | RU |

Hersteller  
Manufacturer  
Fabricant  
Produttore  
Fabricante  
Fabricante  
Производитель

MEDENTiKA® GmbH  
Hammweg 8-10  
76549 Hügelsheim/Germany  
Tel: +49-7229-69912-10  
info@medentika.de  
www.medentika.com

## MRT-Sicherheitsinformationen

### MR Conditional

Bedingt MR-sicher

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## MRI Safety Information

### MR Conditional

Non-clinical testing and MRI simulations were performed to evaluate the dental implant and abutment system offered by MEDENTiKA® GmbH. Non-clinical testing demonstrates that these products are MR Conditional. A patient with an implant or abutment from the MEDENTiKA® Dental System can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5T and 3T, only
- Maximum spatial gradient magnetic field of 4,000 gauss/cm (40 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg and head average SAR of 3.2 W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode

Under the scan conditions defined, the implants and abutments, as well as clinically relevant implant constructs, from the MEDENTiKA® Dental system are expected to produce a maximum temperature rise of 4.9°C after 15-minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the implants or abutments from the MEDENTiKA® Dental System extend approximately 10 mm from this device when imaged with a gradient echo pulse sequence and a 3T MR system.

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